-  (EN) **Molecular Detection Assay 2 - *Listeria monocytogenes***
-  (FR) **Kit de détection moléculaire *Listeria monocytogenes* version 2**
-  (DE) **Molekulare Detektion 2 – *Listeria monocytogenes* Nachweis**
-  (IT) **Analisi molecolare di seconda generazione per il rilevamento di *Listeria monocytogenes***
-  (ES) **Ensayo de Detección Molecular 2 para *Listeria monocytogenes***
-  (NL) **Moleculaire detectie assay 2 *Listeria monocytogenes***
-  (SV) **Molecular Detection Assay 2 – *Listeria monocytogenes***
-  (DA) **Molekylær Detektions Analyse 2 - *Listeria monocytogenes***
-  (NO) **Molekylær deteksjonstest 2 for *Listeria monocytogenes***
-  (FI) **Molekyläärinen testisetti 2 - *Listeria monocytogenes***
-  (PT) **Ensaio para Detecção Molecular de *Listeria monocytogenes* 2**
-  (EL) **Δοκιμασία Μοριακής Ανίχνευσης 2 – *Listeria monocytogenes***
-  (PL) **Molekularny test do wykrywania 2 — *Listeria monocytogenes***
-  (RU) **Молекулярный анализ 2 - *Listeria monocytogenes***
-  (TR) **Moleküler Tayin Testi *Listeria monocytogenes*-2**
-  (JA) **病原菌検出アッセイ2 - リステリア・モノサイトゲネス用**
-  (ZH) **单核细胞增生李斯特菌分子检测分析 2 - 单核细胞增生李斯特菌**
-  (ZH) **分子检测套组 2 - 单增李斯特菌**
-  (TH) **ชุดทดสอบเชื้อระดับโมเลกุล 2 - ลิสทีเรีย โมโนไซโตเจเนส**
-  (KO) **Molecular Detection Assay 2 - 리스테리아 모노사이토제네스**
-  (ID) **Deteksi Molekuler untuk Pengujian *Listeria monocytogenes* 2**

2

Product Instructions

Molecular Detection Assay 2 - *Listeria monocytogenes*

Product Description and Intended Use

3M™ Molecular Detection Assay 2 - *Listeria monocytogenes* is used with the 3M™ Molecular Detection System for the rapid and specific detection of *Listeria monocytogenes* in enriched food and environmental samples. The 3M Molecular Detection Assays use loop-mediated isothermal amplification to rapidly amplify nucleic acid sequences with high specificity and sensitivity, combined with bioluminescence to detect the amplification. Presumptive positive results are reported in real-time while negative results are displayed after the assay is completed. Presumptive positive results should be confirmed using your preferred method or as specified by local regulations^(1, 2, 3).

The 3M Molecular Detection Assay 2 - *Listeria monocytogenes* is intended for use in a laboratory environment by professionals trained in laboratory techniques. 3M has not documented the use of this product in industries other than food or beverage. For example, 3M has not documented this product for testing water, pharmaceutical, cosmetic, clinical or veterinary samples. The 3M Molecular Detection Assay - 2 *Listeria monocytogenes* has not been evaluated with all possible testing protocols or with all possible strains of bacteria.

As with all test methods, the source, formulation and quality of enrichment medium can influence the results. Factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may also influence results. 3M recommends evaluation of the method including enrichment medium, in the user's environment using a sufficient number of samples with particular foods and microbial challenges to ensure that the method meets the user's criteria.

3M has evaluated the 3M Molecular Detection Assay 2 - *Listeria monocytogenes* with Demi-Fraser Broth containing Ferric Ammonium Citrate and Fraser Broth containing Ferric Ammonium Citrate as needed. A typical formulation of this medium follows below.

Demi-Fraser Broth Base Typical Formula	(g/L)	Fraser Broth Base Typical Formula	(g/L)
Sodium Chloride	20 g	Sodium Chloride	20 g
Sodium Phosphate, dibasic, anhydrous*	9.6 g	Sodium Phosphate, dibasic, anhydrous*	9.6 g
Beef Extract	5.0 g	Beef Extract	5.0 g
Pancreatic Digest of Casein	5.0 g	Pancreatic Digest of Casein	5.0 g
Peptic Digest of Animal Tissue	5.0 g	Peptic Digest of Animal Tissue	5.0 g
Yeast Extract	5.0 g	Yeast Extract	5.0 g
Lithium Chloride	3.0 g	Lithium Chloride	3.0 g
Potassium Phosphate, monobasic	1.35 g	Potassium Phosphate, monobasic	1.35 g
Esculin	1.0 g	Esculin	1.0 g
Acriflavin HCl	0.0125 g	Acriflavin HCl	0.025 g
Nalidixic Acid	0.01 g	Nalidixic Acid	0.02 g

* Substitute: Sodium Phosphate, dibasic, dihydrate 12.0 g

Fraser Broth Supplement

(Ingredients per 10 mL vial. One vial is added to one liter of basal medium.)

Ferric Ammonium Citrate 0.5g/10mL

Final pH 7.2 ± 0.2 at 25°C

The 3M™ Molecular Detection Instrument is intended for use with samples that have undergone heat treatment during the assay lysis step, which is designed to destroy organisms present in the sample. Samples that have not been properly heat treated during the assay lysis step may be considered a potential biohazard and should NOT be inserted into the 3M Molecular Detection Instrument.

3M Food Safety is certified to ISO (International Organization for Standardization) 9001 for design and manufacturing.

The 3M Molecular Detection Assay 2 - *Listeria monocytogenes* test kit contains 96 tests, described in Table 1.

Table 1. 3M Molecular Detection Assay Kit Components

Item	Identification	Quantity	Contents	Comments
3M™ Lysis Solution (LS)	Pink solution in clear tubes	96 (12 strips of 8 tubes)	580 µL of LS per tube	Racked and ready to use
3M™ Molecular Detection Assay 2 - <i>Listeria monocytogenes</i> Reagent Tubes	Yellow tubes	96 (12 strips of 8 tubes)	Lyophilized specific amplification and detection mix	Ready to use
Extra caps	Yellow caps	96 (12 strips of 8 caps)		Ready to use
3M™ Reagent Control (RC)	Clear flip-top tubes	16 (2 pouches of 8 individual tubes)	Lyophilized control DNA, amplification and detection mix	Ready to use
Quick Start Guide		1		

The Negative Control, not provided in the kit, is sterile enrichment medium, e.g., Demi-Fraser Broth. Do not use water as a Negative Control.

Safety

The user should read, understand and follow all safety information in the instructions for the 3M Molecular Detection System and the 3M Molecular Detection Assay 2 - *Listeria monocytogenes*. Retain the safety instructions for future reference.

⚠ WARNING: Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

NOTICE: Indicates a potentially hazardous situation which, if not avoided, could result in property damage.

⚠ WARNING

Do not use the 3M Molecular Detection Assay 2 - *Listeria monocytogenes* in the diagnosis of conditions in humans or animals.

The 3M Molecular Detection Assay 2 - *Listeria monocytogenes* method may generate *Listeria monocytogenes* to levels sufficient to cause stillbirths and fatalities in pregnant women and the immunocompromised, if exposed.

The user must train its personnel in current proper testing techniques: for example, Good Laboratory Practices, ISO/IEC 17025⁽⁴⁾, or ISO 7218⁽⁵⁾.

To reduce the risks associated with a false-negative result leading to the release of contaminated product:

- Follow the protocol and perform the tests exactly as stated in the product instructions.
- Store the 3M Molecular Detection Assay 2 - *Listeria monocytogenes* as indicated on the package and in the product instructions.
- Always use the 3M Molecular Detection Assay 2 - *Listeria monocytogenes* by the expiration date.
- Use the 3M Molecular Detection Assay 2 - *Listeria monocytogenes* for food and environmental samples that have been validated internally or by a third party.
- Use the 3M Molecular Detection Assay 2 - *Listeria monocytogenes* only for surfaces, sanitizers, protocols and bacterial strains that have been validated internally or by a third party.
- For an environmental sample containing Neutralizing Buffer with aryl sulfonate complex, perform a 1:2 dilution before testing (1 part sample into 1 part sterile enrichment broth). Another option is to transfer 10 µL of the NB enrichment into the 3M Lysis Solution tubes. 3M™ Sample Collection products which include Neutralizing Buffer with aryl sulfonate complex: BPPFV10NB, RS96010NB, RS9604NB, SSL10NB, SSL10NB2G, HS10NB, HS10NB2G and HS2410NB2G.

To reduce the risks associated with exposure to chemicals and biohazards:

- It is strongly recommended that female laboratory staff be informed of the risk to a developing fetus resulting from infection of the mother through exposure to *Listeria monocytogenes*.
- Perform pathogen testing in a properly equipped laboratory under the control of trained personnel. Incubated enrichment media and equipment or surfaces that have come into contact with incubated enrichment media may contain pathogens at levels sufficient to cause risk to human health.
- Always follow standard laboratory safety practices, including wearing appropriate protective apparel and eye protection while handling reagents and contaminated samples.
- Avoid contact with the contents of the enrichment media and reagent tubes after amplification.
- Dispose of enriched samples according to current industry standards.
- Samples that have not been properly heat treated during the assay lysis step may be considered a potential biohazard and should NOT be inserted into the 3M Molecular Detection Instrument.

To reduce the risks associated with cross-contamination while preparing the assay:

- Always wear gloves (to protect the user and prevent introduction of nucleases).

To reduce the risks associated with environmental contamination:

- Follow current industry standards for disposal of contaminated waste.

To reduce the risks associated with exposure to hot liquids:

- Do not exceed the recommended temperature setting on heater.
- Do not exceed the recommended heating time.
- Use an appropriate, calibrated thermometer to verify the 3M™ Molecular Detection Heat Block Insert temperature (e.g., a partial immersion thermometer or digital thermocouple thermometer, not a total immersion thermometer.) The thermometer must be placed in the designated location in the 3M Molecular Detection Heat Block Insert.

NOTICE

To reduce the risks associated with cross-contamination while preparing the assay:

- Change gloves prior to reagent pellet hydration.
- Use of sterile, aerosol barrier (filtered), molecular biology grade pipette tips is recommended.
- Use a new pipette tip for each sample transfer.
- Use Good Laboratory Practices to transfer the sample from the enrichment to the lysis tube. To avoid pipettor contamination, the user may choose to add an intermediate transfer step. For example, the user can transfer each enriched sample into a sterile tube.
- Use a molecular biology workstation containing germicidal lamp where available. Periodically decontaminate laboratory benches and equipment (pipettes, cap/decap tools, etc.) with a 1- 5% (v:v in water) household bleach solution or DNA removal solution.

To reduce the risks associated with a false-positive result:

- Never open reagent tubes post amplification.
- Always dispose of the contaminated tubes by soaking in a 1-5% (v:v in water) household bleach solution or equivalent solution for 1 hour and away from the assay preparation area.
- Never autoclave reagent tubes post amplification.

Consult the Safety Data Sheet for additional information and local regulations for disposal.

If you have questions about specific applications or procedures, please visit our website at www.3M.com/foodsafety or contact your local 3M representative or distributor.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.3M.com/foodsafety, or contact your local 3M representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user’s responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user’s criteria.

It is also the user’s responsibility to determine that any test methods and results meet its customers’ and suppliers’ requirements.

As with any test method, results obtained from use of any 3M Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

To help customers evaluate the method for various food matrices, 3M has developed the 3M™ Molecular Detection Matrix Control kit. When needed, use the Matrix Control (MC) to determine if the matrix has the ability to impact the 3M Molecular Detection Assay 2 - *Listeria monocytogenes* results. Test several samples representative of the matrix, i.e., samples obtained from different origin, during any validation period when adopting the 3M method or when testing new or unknown matrices or matrices that have undergone raw material or process changes.

A matrix can be defined as a type of product with intrinsic properties such as composition and process. Differences between matrices may be as simple as the effects caused by differences in their processing or presentation for example, raw vs. pasteurized; fresh vs. dried, etc.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, 3M DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any 3M Food Safety Product is defective, 3M or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify 3M within sixty days of discovery of any suspected defects in a product and return it to 3M. Please call Customer Service (1-800-328-1671 in the U.S.) or your official 3M Food Safety representative for a Returned Goods Authorization.

Limitation of 3M Liability

3M WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall 3M's liability under any legal theory exceed the purchase price of the product alleged to be defective.

Storage and Disposal

Store the 3M Molecular Detection Assay 2 - *Listeria monocytogenes* at 2-8°C. Do not freeze. Keep kit away from light during storage. After opening the kit, check that the foil pouch is undamaged. If the pouch is damaged, do not use. After opening, unused reagent tubes should always be stored in the re-sealable pouch with the desiccant inside to maintain stability of the lyophilized reagents. Store resealed pouches at 2-8°C for no longer than 60 days.

Do not use 3M Molecular Detection Assay 2 - *Listeria monocytogenes* past the expiration date. Expiration date and lot number are noted on the outside label of the box. After use, the enrichment medium and the 3M Molecular Detection Assay 2 - *Listeria monocytogenes* tubes can potentially contain pathogenic materials. When testing is complete, follow current industry standards for the disposal of contaminated waste. Consult the Safety Data Sheet for additional information and local regulations for disposal.

Instructions for Use

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Periodically decontaminate laboratory benches and equipment (pipettes, cap/decap tools, etc.) with a 1-5% (v:v in water) household bleach solution or DNA removal solution.

The user should complete the 3M Molecular Detection System operator qualification (OQ) training, as described in the "Installation Qualification (IQ) / Operational Qualification (OQ) Protocols and Instructions for 3M Molecular Detection System" document⁽⁶⁾.

See Section "Specific Instructions for validated methods" for specific requirements:

Table 3 for enrichment protocols according to AOAC® *Official Method of Analysis*SM **2016.08** and AOAC® *Performance Tested Method*SM Certificate **#081501**.

Table 4 for enrichment protocols according to NF Validation certificate **3M 01/15-09/16**.

Sample Enrichment

Tables 2, 3 or 4 present guidance for the enrichment of food and environmental samples. It is the user's responsibility to validate alternate sampling protocols or dilution ratios to ensure this test method meets the user's criteria.

Foods

1. Allow the Demi-Fraser Broth enrichment medium (includes ferric ammonium citrate) to equilibrate to ambient laboratory temperature.
2. Aseptically combine the enrichment medium and sample according to Tables 2, 3 or 4. For all meat and highly particulate samples, the use of filter bags is recommended.
3. Homogenize thoroughly by blending, stomaching, or hand mixing for 2 ± 0.2 minutes. Incubate at $37 \pm 1^\circ\text{C}$ according to Tables 2, 3 or 4.
4. For raw dairy products (see Tables 2, 3 or 4), transfer 0.1 mL of the primary enrichment into 10 mL of Fraser Broth. Incubate at $37 \pm 1^\circ\text{C}$ for 20-24 hours.

Environmental samples

Sample collection devices can be a sponge hydrated with a neutralizing solution to inactivate the effects of the sanitizers. 3M recommends the use of a biocide-free cellulose sponge. Neutralizing solution can be Dey-Engley (D/E) Neutralizing Broth or Letheen Broth. It is recommended to sanitize the area after sampling.

WARNING: Should you select to use Neutralizing Buffer (NB) that contains aryl sulfonate complex as the hydrating solution for the sponge, it is required to perform a 1:2 dilution (1 part sample into 1 part sterile enrichment broth) of the enriched environmental sample before testing in order to reduce the risks associated with a false-negative result leading to the release of contaminated product. Another option is to transfer 10 μL of the neutralizing buffer enrichment into the 3M Lysis Solution tubes.

The recommended size of the sampling area to verify the presence or absence of the pathogen on the surface is at least 100 cm² (10 cm x 10 cm or 4" x 4"). When sampling with a sponge, cover the entire area going in two directions (left to right then up and down) or collect environmental samples following your current sampling protocol or according to the FDA BAM⁽¹⁾, USDA FSIS MLG⁽²⁾ or ISO 18593⁽⁷⁾ guidelines.

1. Allow the Demi-Fraser Broth enrichment medium (includes ferric ammonium citrate) to equilibrate to ambient laboratory temperature.

2. Aseptically combine the enrichment medium and sample according to Tables 2, 3 or 4.
3. Homogenize thoroughly by blending, stomaching, or hand mixing for 2 ± 0.2 minutes. Incubate at $37 \pm 1^\circ\text{C}$ for 24-30 hours according to Tables 2, 3, or 4.

Table 2: General enrichment protocols at $37 \pm 1^\circ\text{C}$ using Demi-Fraser Broth and Fraser Broth as needed.

Sample Matrix	Sample Size	Enrichment Broth Volume (mL)	Enrichment Time (hours)				
Heat-processed, cooked, cured meats, poultry, seafood and fish	25 g	225	24-30				
Heat-processed / pasteurized dairy products							
Produce and vegetables							
Multi-component foods							
Environmental samples	1 sponge	100 or 225	24-30				
	1 swab	10	24-30				
Raw meat, poultry, seafood, fish	25 g	475	28-32				
Sample Matrix	Primary Enrichment (Demi-Fraser Broth) ^(a)			Secondary Enrichment (Fraser Broth) ^(a)			Sample Analysis Volume ^(b)
	Sample Size	Enrichment Broth Volume (mL)	Enrichment Time (hours)	Sample Size	Enrichment Temperature ($^\circ\text{C}$)	Enrichment Time (hours)	
Raw dairy products	25 g	225	20-24	Transfer 0.1 mL into 10 mL Fraser Broth	37 ± 1	20-24	10 μL

- (a) Demi-Fraser and Fraser Broth should always be supplemented with Fraser Broth supplement (ferric ammonium citrate) during primary or secondary enrichment.
- (b) Volume of sample transferred to the 3M Lysis Solution tubes. Refer to Step 4.6 of Lysis section.

Specific Instructions for Validated Methods
AOAC® Official Method of AnalysisSM 2016.08
AOAC® Performance Tested MethodSM #081501

In AOAC Research Institute OMASM and PTMSM programs, the 3M Molecular Detection Assay 2 - *Listeria monocytogenes* was found to be an effective method for the detection of *Listeria monocytogenes*. The matrices tested in the study are shown in Table 3.

Table 3. Enrichment protocols using Demi-Fraser Broth^(a) at 37 ± 1 °C according to AOAC OMASM 2016.08 and PTMSM Certificate #081501.

Sample Matrix		Sample Size	Enrichment Broth Volume (mL)	Enrichment Time (hours)
Beef hot dogs, queso fresco, vanilla ice cream, 4% milk fat cottage cheese, 3% chocolate whole milk, romaine lettuce, bagged raw spinach, cold smoked salmon		25 g	225	24-30
Raw chicken		25 g	475	28-32
Deli turkey		125 g	1125	24-30
Cantaloupe ^(b)		Whole melon	Enough volume to allow melon to float	26-30
Environmental samples:	Stainless steel	1 sponge	225	24-30
	Sealed concrete	1 sponge	100	24-30
	Plastic ^(c)	1 swab	10	24-30

All samples for the AOAC validation were homogenized by stomaching unless otherwise noted.

- (a) Demi-Fraser and Fraser Broth should always be supplemented with Fraser Broth supplement (ferric ammonium citrate) during primary or secondary enrichment.
- (b) Homogenize sample by hand mixing.
- (c) Homogenize sample by vortexing.

NF VALIDATION by AFNOR Certification

3M 01/15-09/16

ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS

<http://nf-validation.afnor.org/en>

For more information about end of validity, please refer to NF VALIDATION certificate available on the website mentioned above.

NF VALIDATION certified method in compliance with ISO 16140-2⁽⁸⁾ in comparison to ISO 11290⁽³⁾

Scope of the validation: All human food and environmental samples (excluding primary production samples)

Sample preparation: Samples should be prepared according to EN ISO 11290-1⁽³⁾ and EN ISO 6887⁽⁹⁾

Software version: See certificate

Table 4: Enrichment protocols according to NF VALIDATION certified method 3M 01/15-09/16 at 37 ± 1°C using Demi-Fraser Broth and Fraser Broth as needed.

General Protocol	Sample Size	Enrichment Broth Volume (mL)	Enrichment Temperature (±1°C)	Enrichment Time (hours)	Sample Analysis Volume ^(a)	Accepted Interruption Point
All food samples (excepted raw meats, raw seafood, and raw dairy products)	25 g	225	37	24-30	20 µL	<ul style="list-style-type: none"> • Demi-Fraser Broth up to 72 hours • Lysate at -20°C • Lysate at 4°C up to 72 hours
Environmental samples	25 g, 1 swab, or 1 wipe					

Specific Protocol	Primary Enrichment (Demi-Fraser Broth) ^(b)				Secondary Enrichment (Fraser Broth) ^(b)				Accepted Interruption Point
	Sample Size	Enrichment Broth Volume (mL)	Enrichment Temperature ($\pm 1^{\circ}\text{C}$)	Enrichment Time (hours)	Sample Size	Enrichment Temperature ($\pm 1^{\circ}\text{C}$)	Enrichment Time (hours)	Sample Analysis Volume ^(c)	
Raw meats, raw seafood, and raw dairy products	25 g	225	37	20-24	Transfer 0.1 mL into 10 mL Fraser Broth	37	20-24	10 μL	<ul style="list-style-type: none"> Fraser Broth up to 72 hours Lysate at -20°C Lysate at 4°C up to 72 hours

(a) Volume of sample transferred to the 3M Lysis Solution tubes. Refer to Step 4.6 of Lysis section.

(b) Demi-Fraser and Fraser Broth should always be supplemented with Fraser Broth supplement (ferric ammonium citrate) during primary or secondary enrichment.

(c) Volume of sample transferred to the 3M Lysis Solution tubes. Refer to Step 4.6 of Lysis section.

NOTE: Samples larger than 25 g have not been tested in the NF validation study.

Preparation of the 3M™ Molecular Detection Speed Loader Tray

1. Wet a cloth with a 1-5% (v:v in water) household bleach solution or equivalent solution and wipe the 3M™ Molecular Detection Speed Loader Tray.
2. Rinse the 3M Molecular Detection Speed Loader Tray with water.
3. Use a disposable towel to wipe the 3M Molecular Detection Speed Loader Tray dry.
4. Ensure the 3M Molecular Detection Speed Loader Tray is dry before use.

Preparation of the 3M™ Molecular Detection Chill Block Insert

Place the 3M™ Molecular Detection Chill Block directly on the laboratory bench; Use the chill block at ambient laboratory temperature ($20\text{-}25^{\circ}\text{C}$).

Preparation of the 3M™ Molecular Detection Heat Block Insert

Place the 3M™ Molecular Detection Heat Block Insert in a dry double block heater unit. Turn on the dry block heater unit and set the temperature to allow the 3M Molecular Detection Heat Block Insert to reach and maintain a temperature of $100 \pm 1^{\circ}\text{C}$.

NOTE: Depending on the heater unit, allow approximately 30 minutes for the 3M Molecular Detection Heat Block Insert to reach temperature. Using an appropriate, calibrated thermometer (e.g., a partial immersion thermometer or digital thermocouple thermometer, **not** a total immersion thermometer) placed in the designated location, verify that the 3M Molecular Detection Heat Block Insert is at $100 \pm 1^{\circ}\text{C}$.


Preparation of the 3M™ Molecular Detection Instrument

1. Launch the 3M™ Molecular Detection System Software and log in. Contact your 3M Food Safety representative to ensure you have the most updated version of the software.
2. Turn on the 3M Molecular Detection Instrument.
3. Create or edit a run with data for each sample. Refer to the 3M Molecular Detection System User Manual for details.

NOTE: The 3M Molecular Detection Instrument must reach Ready State before inserting the 3M Molecular Detection Speed Loader Tray with reaction tubes. This heating step takes approximately 20 minutes and is indicated by an ORANGE light on the instrument's status bar. When the instrument is ready to start a run, the status bar will turn GREEN.


Lysis

Remove the bottom of 3M Lysis Solution Rack with a screwdriver or spatula before placing in in the 3M Molecular Detection Heat Block Insert.


1. Allow the lysis solution (3M Lysis Solution) tubes to warm up by setting the rack at room temperature (20-25°C) overnight (16-18 hours). Alternatives to equilibrate the 3M Lysis Solution tubes to room temperature are to set the 3M Lysis Solution tubes on the laboratory bench for at least 2 hours, incubate the 3M Lysis Solution tubes in a 37 ± 1°C incubator for 1 hour or place them in a dry double block heater for 30 seconds at 100°C.
2. Invert the capped tubes to mix. Proceed to next step within 4 hrs.
3. Remove the enrichment broth from the incubator.
4. One 3M Lysis Solution tube is required for each sample and the Negative Control (NC) sample (sterile enrichment medium).
 - 4.1 3M Lysis Solution tube strips can be cut to desired 3M Lysis Solution tube number. Select the number of individual 3M Lysis Solution tubes or 8-tube strips needed. Place the 3M Lysis Solution tubes in an empty rack.
 - 4.2 To avoid cross-contamination, decap one 3M Lysis Solution tubes strip at a time and use a new pipette tip for each transfer step.
 - 4.3 Transfer enriched samples to 3M Lysis Solution Tubes as described below:

Transfer each enriched sample into an individual 3M Lysis Solution tube **first**. Transfer the NC **last**.
 - 4.4 Use the 3M™ Molecular Detection Cap/Decap Tool-Lysis to decap one 3M Lysis Solution tube strip – one strip at a time.
 - 4.5 Discard the 3M Lysis Solution tube cap – if lysate will be retained for retest, place the caps into a clean container for re-application after lysis.
 - 4.5.1 For processing of retained lysate, see Appendix A.
 - 4.6 Transfer 20 µL of sample into a 3M Lysis Solution tube unless otherwise indicated in Protocols from Tables 2, 3, or 4 (e.g., raw dairy products or when using environmental samples with neutralizing buffer use 10 µL).

5. Repeat step 4.4 to 4.6 until each individual sample has been added to a corresponding 3M Lysis Solution tube in the strip.
6. When all samples have been transferred, transfer 20 µL of NC (sterile enrichment medium, e.g., Demi-Fraser Broth) into a 3M Lysis Solution tube. Do not use water as a NC.
7. Verify that the temperature of the 3M Molecular Detection Heat Block Insert is at 100 ± 1°C.
8. Place the uncovered rack of 3M Lysis Solution tubes in the 3M Molecular Detection Heat Block Insert and heat for 15 ± 1 minutes. During heating, the 3M Lysis Solution will change from pink (cool) to yellow (hot).
 - 8.1 Samples that have not been properly heat treated during the assay lysis step may be considered a potential biohazard and should NOT be inserted into the 3M Molecular Detection Instrument.
9. Remove the uncovered rack of 3M Lysis Solution tubes from the heating block and allow to cool in the 3M Molecular Detection Chill Block Insert at least 5 minutes and a maximum of 10 minutes. The 3M Molecular Chill Block Insert used at ambient temperature should sit directly on the laboratory bench. When cool, the lysis solution will revert to a pink color.

10. Remove the rack of 3M Lysis Solution tubes from the 3M Molecular Detection Chill Block Insert.


Amplification

1. One 3M Molecular Detection Assay 2 - *Listeria monocytogenes* Reagent Tube is required for each sample and the NC.
 - 1.1 Reagent tube strips can be cut to desired tube number. Select the number of individual 3M Molecular Detection Assay 2 - *Listeria monocytogenes* Reagent Tubes or 8-tube strips needed.
 - 1.2 Place 3M Molecular Detection Assay 2 - *Listeria monocytogenes* Reagent Tubes in an empty rack.
 - 1.3 Avoid disturbing the reagent pellets from the bottom of the tubes.
2. Select 1 3M Reagent Control tube and place in rack.
3. To avoid cross-contamination, decap one 3M Molecular Detection Assay 2 - *Listeria monocytogenes* Reagent Tube strip at a time and use a new pipette tip for each transfer step.
4. Transfer lysate to 3M Molecular Detection Assay 2 - *Listeria monocytogenes* Reagent Tubes and 3M Reagent Control tube as described below:

Transfer each sample lysate into individual 3M Molecular Detection Assay 2 - *Listeria monocytogenes* Reagent Tubes **first** followed by the NC. Hydrate the 3M Reagent Control tube **last**.

5. Use the 3M™ Molecular Detection Cap/Decap Tool-Reagent to decap the 3M Molecular Detection Assay 2 - *Listeria monocytogenes* Reagent Tubes –one Reagent tube strip at a time. Discard cap.
 - 5.1 **Transfer 20 µL of Sample lysate from the upper ½ of the liquid (avoid precipitate) in the 3M Lysis Solution tube into corresponding 3M Molecular Detection Assay 2 - *Listeria monocytogenes* Reagent Tube. Dispense at an angle to avoid disturbing the pellets. Mix by gently pipetting up and down 5 times.**
 - 5.2 Repeat step 5.1 until individual Sample lysate has been added to a corresponding 3M Molecular Detection Assay 2 - *Listeria monocytogenes* Reagent Tube in the strip.
 - 5.3 Cover the 3M Molecular Detection Assay 2 - *Listeria monocytogenes* Reagent Tubes with the provided extra caps and use the rounded side of the 3M Molecular Detection Cap/Decap Tool-Reagent to apply pressure in a back and forth motion ensuring that the cap is tightly applied.
 - 5.4 Repeat step 5.1 to 5.3 as needed, for the number of samples to be tested.
 - 5.5 When all sample lysates have been transferred, repeat steps 5.1 to transfer 20 µL of NC lysate into a 3M Molecular Detection Assay 2 - *Listeria monocytogenes* Reagent Tube.
 - 5.6 Transfer **20 µL of NC lysate into a 3M Reagent Control tube**. Dispense at an angle to avoid disturbing the pellets. Mix by gently pipetting up and down 5 times.
6. Load capped tubes into a clean and decontaminated 3M Molecular Detection Speed Loader Tray. Close and latch the 3M Molecular Detection Speed Loader Tray lid.

7. Review and confirm the configured run in the 3M Molecular Detection System Software.
8. Click the Start button in the software and select instrument for use. The selected instrument's lid automatically opens.

9. Place the 3M Molecular Detection Speed Loader Tray into the 3M Molecular Detection Instrument and close the lid to start the assay. Results are provided within 75 minutes, although positives may be detected sooner.
10. After the assay is complete, remove the 3M Molecular Detection Speed Loader Tray from the 3M Molecular Detection Instrument and dispose of the tubes by soaking in a 1-5% (v:v in water) household bleach solution or equivalent solution for 1 hour and away from the assay preparation area.

NOTICE: To minimize the risk of false positives due to cross-contamination, never open reagent tubes containing amplified DNA. This includes 3M Molecular Detection Assay 2 - *Listeria monocytogenes* Reagent Tubes, 3M Reagent Control and Matrix Control tubes. Always dispose of sealed reagent tubes by soaking in a 1-5% (v:v in water) household bleach solution or equivalent solution for 1 hour and away from the assay preparation area.

Results and Interpretation

An algorithm interprets the light output curve resulting from the detection of the nucleic acid amplification. Results are analyzed automatically by the software and are color-coded based on the result. A Positive or Negative result is determined by analysis of a number of unique curve parameters. Presumptive positive results are reported in real-time while Negative and Inspect results will be displayed after the run is completed.

Presumptive positive samples should be confirmed as per the laboratory standard operating procedures or by following the appropriate reference method confirmation^(1,2,3), beginning with transfer from the primary enrichment to secondary enrichment broth (if applicable), followed by subsequent plating and confirmation of isolates using appropriate biochemical and serological methods.

NOTE: Even a negative sample will not give a zero reading as the system and 3M Molecular Detection Assay 2 - *Listeria monocytogenes* amplification reagents have a “background” relative light unit (RLU) reading.

In the rare event of any unusual light output, the algorithm labels this as “Inspect.” 3M recommends the user to repeat the assay for any Inspect samples. If the result continues to be Inspect, proceed to confirmation test using your preferred method or as specified by local regulations.

Table 5: Result symbols for samples and interpretation as provided by the 3M Molecular Detection System Software.

Well Type	Well Result Symbol	Result	Interpretation
Sample		Positive	The sample is presumptive positive for the target pathogen.
Sample		Negative	The sample is negative for the target pathogen.
Sample		Inhibited	The sample matrix was inhibitory to the assay. A re-test may be required. Refer to the troubleshooting section and the assay kit Product Instructions for more information.
Sample		Inspect	The presence or absence of the target pathogen was indeterminate. A re-test may be required. Refer to the troubleshooting section and the assay kit Product Instructions for more information.
Sample		Error	No bioluminescence was detected. A re-test may be required. Refer to the troubleshooting section and the assay kit Product Instructions for more information.

Confirmation of Results According to the NF VALIDATION Certified Method

Option 1: Using the ISO 11290⁽³⁾ standard starting from Demi-Fraser enrichment

Option 2: By isolation on PALCAM agar or on chromogenic agar forming part of a NF VALIDATION certified method for detection of the *Listeria monocytogenes*. The presence of characteristic colonies is sufficient to confirm the presence of *Listeria monocytogenes*.

Option 3: Using nucleic acid probes as described in EN ISO 7218⁽⁵⁾ standard, performed on isolated colonies, from selective agar (see Options 1 or 2).

Option 4: Using any other method certified NF VALIDATION, the principle of which must be different from 3M Molecular Detection Assay 2 - *Listeria monocytogenes*. The complete protocol described for this second validated method must be used. All steps prior to the start of confirmation must be common to both methods.

In the event of discordant results (presumptive positive with the alternative method, non-confirmed by one of the means described above) the laboratory must follow the necessary steps to ensure the validity of the result obtained.

In the event of discordant results (presumptive positive with the 3M Molecular Detection Assay 2 - *Listeria monocytogenes*, non-confirmed by one of the means described above, the laboratory must follow the necessary steps to ensure the validity of the results obtained.

If you have questions about specific applications or procedures, please visit our website at www.3M.com/foodsafety or contact your local 3M representative or distributor.

Appendix A. Protocol Interruption: Storage and re-testing of heat-treated lysates.

1. To store a heat-treated lysate, re-cap the lysis tube with a clean cap (see “Lysis”, 4.5)
2. Store at 4 to 8°C for up to 72 hours.
3. Prepare a stored sample for amplification by inverting 2-3 times to mix.
4. Decap the tubes.
5. Place the mixed lysate tubes on 3M Molecular Detection Heat Block Insert and heat at $100 \pm 1^\circ\text{C}$ for 5 ± 1 minutes.
6. Remove the rack of 3M Lysis Solution tubes from the heating block and allow to cool in the 3M Molecular Detection Chill Block Insert at least 5 minutes and a maximum of 10 minutes.
7. Continue the protocol at the ‘Amplification’ section detailed above.

References:

1. US Food and Drug Administration Bacteriological Analytical Manual. Chapter 10: Detection of *Listeria monocytogenes* in foods and environmental samples, and enumeration of *Listeria monocytogenes* in foods. October 31, 2017.
2. US Department of Agriculture (USDA) FSIS Microbiology Laboratory Guidebook 8.11. Isolation and identification of *Listeria monocytogenes* from red meat, poultry, ready-to-eat siluriformes (fish) and egg products, and environmental samples. January 2, 2019.
3. NF EN ISO 11290-1. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp. – Part 1: Detection method.
4. NF EN ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. NF EN ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. 3M Installation Qualification (IQ) / Operational Qualification (OQ) Protocols and Instructions for 3M Molecular Detection System.
7. NF EN ISO 18593. Microbiology of the food chain – Horizontal methods for surface sampling.
8. NF EN ISO 16140-2 Microbiology of the food chain – Method validation – Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method.
9. NF EN ISO 6887. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.

Explanation of Product Label Symbols

www.3M.com/foodsafety/symbols

3M Food Safety

3M United States

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-800-328-6553

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-563-2921

3M Latin America

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-954-340-8263

3M Europe and MEA

3M Deutschland GmbH
Carl-Schurz - Strasse 1
D41453 Neuss/Germany
+49-2131-14-3000

3M United Kingdom PLC

Morley Street,
Loughborough
Leicestershire
LE11 1EP
United Kingdom
+(44) 1509 611 611

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3M Australia

Bldg A, 1 Rivett Road
North Ryde, NSW 2113
Australia
61 1300 363 878

3M Company
2510 Conway Ave
St. Paul, MN 55144 USA
www.3M.com/foodsafety

3M is a trademark of 3M.
Used under license in Canada. © 2020, 3M.
Unauthorized use prohibited. All rights reserved.

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Demi Fraser Broth Base 성분표	(g/L)	Fraser Broth Base의 성분표	(g/L)
Sodium Chloride	20g	Sodium Chloride	20g
Sodium Phosphate, dibasic, anhydrous*	9.6g	Sodium Phosphate, dibasic, anhydrous*	9.6g
Beef Extract	5.0g	Beef Extract	5.0g
Pancreatic Digest of Casein	5.0g	Pancreatic Digest of Casein	5.0g
Peptic Digest of Animal Tissue	5.0g	Peptic Digest of Animal Tissue	5.0g
Yeast Extract	5.0g	Yeast Extract	5.0g
Lithium Chloride	3.0g	Lithium Chloride	3.0g
Potassium Phosphate, monobasic	1.35g	Potassium Phosphate, monobasic	1.35g
Esculin	1.0g	Esculin	1.0g
Acridine HCl	0.0125g	Acridine HCl	0.025g
Nalidixic Acid	0.01g	Nalidixic Acid	0.02g

* 대체: Sodium Phosphate, dibasic, dihydrate 12.0g

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Referensi:

1. US Food and Drug Administration Bacteriological Analysis Manual. Chapter 10: Detection of *Listeria monocytogenes* in foods and environmental samples, and enumeration of *Listeria monocytogenes* in foods. October 31, 2017.
2. US Department of Agriculture (USDA) FSIS Microbiology Laboratory Guidebook 8.11. Isolation and identification of *Listeria monocytogenes* from red meat, poultry, ready-to-eat siluriformes (fish) and egg products, and environmental samples. January 2, 2019.
3. NF EN ISO 11290-1. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp. – Part 1: Detection method.
4. NF EN ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. NF EN ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. 3M Installation Qualification (IQ) / Operational Qualification (OQ) Protocols and Instructions for 3M Molecular Detection System.
7. NF EN ISO 18593. Microbiology of the food chain – Horizontal methods for surface sampling.
8. NF EN ISO 16140-2 Microbiology of the food chain – Method validation – Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method.
9. NF EN ISO 6887. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.

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3M Health Care Limited